

SPURTEX®

DISINFECTION OF NANOFIBER-BASED PRODUCTS

SPURTEX® disinfection is a liquid disinfectant with a virucidal effect intended for the general public for single and repeated disinfection of nanofiber-based face masks and respirators.

- For disinfecting SPURTEX® nanofiber-based face masks and respirators
- Also suitable as hand disinfection
- Effective and immediate inactivation of viruses and bacteria
- Allows reuse of SPURTEX® nanofiber-based face masks and respirators up to 5 times
- It saves money and the environment
- Odour-free



spray **2–3 times**
from each side



let dry for
24 hours



the process can
be repeated
up to 5 times



✓
IMMEDIATE EFFECT



Czech-made

Instructions for use

1. Wash your hands thoroughly and disinfect them with SPURTEX® disinfection before and immediately after disinfecting the face mask/respirator.
2. Spray all parts of the face mask/respirator evenly and over the entire surface with disinfectant. Apply 2-3 sprays on both sides from a distance of about 15 cm.
3. Then allow the face mask/respirator to dry and air it at room temperature, ideally for 24 hours, but at least overnight.
4. After thorough drying, the face mask/respirator can be used again.
5. The process can be repeated up to 5 times.

Additional information:

The face mask/respirator should be disinfected after each use (if possible to do so). Shake the disinfectant bottle well before use. To reduce the risk of inhaling residual alcohol fumes, do not spray on the face mask/respirator you are currently wearing. It is not recommended to use a dryer for drying, as this could mechanically damage the filtration layer. Increased thermal stress on the face mask/respirator due to drying in a microwave or with a hair dryer is not recommended either. Various means such as ethanol, isopropanol can be used to disinfect the face mask/respirator. The World Health Organization (WHO) recommends the use of alcohol at a concentration of 70 to 90% for disinfection. SPUR a.s. company offers UNISANOL alcohol-based disinfectants. UNISANOL disinfectant affects only the aesthetic form (appearance); it does not affect the functionality of the filter material.

Supplier:
SPUR a.s. | Tomáše Bati 299, Louky, 763 02 Zlín,
Czech Republic
Phone: +420 577 601 112 | e-mail: spurtex@spur.cz |
www.spur.cz | eshop.spur.cz



Producer:
Tomas Bata University in Zlín
(Unisanol recipe – reg. no. MZDR 17099/2020/OBP).

